

Forward and inverse solutions for EEG/MEG source reconstruction

Amit Jaiswal, DSc (Tech.)

Senior Clinical Scientist & Application Lead

MEGIN Oy, Espoo, Finland



July 25, 2025

IIT Mandi, Mandi, HP, India




MEG/EEG signals measurement

MEG and EEG track electric brain activity by measuring the electromagnetic fields generated by neurons

EEG = measuring the potential differences on the scalp

MEG = measuring neuromagnetic fields outside of the head




MEG + EEG




Sensor-level vs. source-level analysis


- Often simple and fast and generally involve fewer steps.
- Doesn't require brain anatomical information.
- Robust but lower spatial accuracy.






- Involves multiple steps, turning the analysis workflow more complex and time-consuming.
- Often require anatomical information of patient's head.
- Various influencing factors.



A typical workflow for EEG source imaging



Forward and inverse models: simplified



Forward model

Inverse model

Need of forward solution for EEG analysis


Forward model/solution is a key ingredient in EEG (or MEG) source localization, enabling *inverse methods* to estimate brain activity from recorded EEG signals.

Purpose: To predict the electrical potentials at scalp electrodes for a known brain source activity, e.g. a dipole.

Inputs: Requires a head model, a source model, and sensors position. Structural MRIs are often utilized for head and source model.

Computation: Solves Maxwell's equations (quasi-static approximation) to estimate how currents propagate through brain tissues to the scalp.

Output: Produces a lead field matrix that maps each source's contribution to each EEG electrode.



Forward model formulation

The EEG forward model can be written as:

$$\mathbf{V} = \mathbf{L} \cdot \mathbf{J}$$

where,

V: Sensor data (EEG signals)


shape: [#channels x times]

L: Leadfield matrix (describes how sources project to sensors)

shape: [#channels × 3 x #sources]


J: Source currents (dipole amplitudes at each source location)

shape: [3 x nsources]




Use of structural MRI from subject's head

- **Need:** To define head compartments and source space.
- T1-weighted structural MRIs are often used.
- Various MRI processing software can be used to segment the MRI and reconstruct the head surface, e.g., FSL, Freesurfer, and SPM.
- Improve localization accuracy; however it poses additional cost.
- Also intersecting surfaces are a major challenge in EEG source imaging.




Volume conductor model *aka.* head model


- **Physical representation of head tissues:** Models the head as a set of compartments with distinct-conductivities.
- **Includes geometry and conductivity:** It accounts for the shape and size of head structures, as well as their anisotropic or isotropic conductivity properties.
- **Realistic or simplified:** Can be modeled with a single or multiple compartments. It may be a simple spherical model and realistic models BEM or FEM.



Single or multiple compartment conductor model?






- Often used and sufficient for MEG




- EEG source reconstruction often requires a multiple compartment head model

Source models

- Dipolar sources, placed within the brain boundary.
- They represent the locations and orientations of the neural sources).
- Could be either a 3D grid (*volume source model*) or constrained to the cortical surface (*cortically-constrained source model*).



Digitization and coregistration



Jaiswal, et al. (2023). On electromagnetic head digitization in MEG and EEG. *Scientific Reports*, 13(1), 3801.



Challenges with using individual MRIs with EEG

- Additional time, cost and logistics.
- Additional data processing.
- Pediatric and metal-implanted patients often excluded
- Occasionally insufficient quality
 - Deformed scalp surface may affect coregistration.
 - Automatic segmentation leads to incorrect / intersecting surfaces





Alternative approaches to individual MRIs

- Sphere model
- Template MRI
- Scaled template MRI



Inner skull surface from individual MRI


Using MRI templates and pseudo-MRI



Jaiswal, et al. (2025). **Pseudo-MRI Engine** for MRI-Free Electromagnetic Source Imaging. *Human Brain Mapping*, 46(2), e70148.


EEG source reconstruction with pseudo-MRI *

- Simulated spikes




	<i>real MRI</i>	<i>pseudo-MRI</i>
Location	43.0, -23.0, 29.0	37.0, -23.0, 24.0
Orientation	-0.4, 0.3, 0.84	-0.4, -0.4, -0.8
Amplitude	89 nAm	91 nAm
GOF	63.3	63.8
Loc. error	15.2 mm	13.1 mm

- Binaural auditory (tone) stimulus



	<i>real MRI</i>	<i>pseudo-MRI</i>
Location	33.0, -2.0, 57.0	27.0, 0.0, 55.0
Orientation	0.1, -0.6, -0.8	0.2, -0.6, -0.8
Amplitude	23 nAm	23 nAm
GOF	67.8	70.0




* Study in progress (Pseudo-MRI in EEG source imaging)

MEGIN

Inverse methods

The inverse methods for computing a source reconstruction can be divided into three categories:

- **Dipole fitting:** using an overdetermined model with a few sources
 - Single dipole fitting
 - Sequential dipole fitting
- **Scanning:** computed independently on each point of a source model
 - dynamic imaging of coherent sources (DICS)
 - linear constrained minimum variance (LCMV)
- **Distributed source modeling:** using an underdetermined distributed source model
 - minimum norm estimation with and without noise regularisation (MNE)
 - minimum norm estimation using eLORETA




Inverse modeling

- **Goal:** To estimate the location, orientation, and strength of neural sources in the brain that produce the recorded EEG signals.
- **Ill-posed problem:** The inverse problem is ill-posed, means there is no unique solution without additional constraints.

Why? Because there are **more unknowns** (sources) than equations (EEG sensors); which make the problem **ill-posed**, as many different J can explain the same V in $V = L.J$

- **Prior constraint:** Prior knowledge (e.g., sources only on cortex, fixed orientations) helps reduce ambiguity and improve accuracy.
- **Different inverse methods:** Various methods including beamforming, sLORETA, eLORETA, MNE, and sparse solvers, follow unique assumptions.





Let's see how we do it in
FieldTrip!

Read MRI data and sensor position

```
mri = ft_read_mri(mri_fname);  
ft_sourceplot([], mri)
```



```
ft_plot_sens(raw_clean.elec, 'label', 'off',...  
            'elecsz', 1, 'elecshape', 'circle',...  
            'facecolor', 'red');
```

```
hsps = ft_read_headshape(eeg_filename);  
ft_plot_headshape(hsps, 'vertexcolor', 'gray',...  
                  'vertexsize', 15)
```




Coregistration


- Interactive mode
- Headshape-base




Segmentation of MRI



brain



skull



scalp


Head model preparation

`% Prepare triangular mesh`




```
cfg      = [];  
cfg.method = 'projectmesh';  
cfg.numvertices = [3000 2000 1000];  
cfg.tissue = {'brain', 'skull', 'scalp'};  
mesh      = ft_prepare_mesh(cfg, segmri);
```

`% compute the 3-compartment conductor model`

```
cfg      = [];  
cfg.method = 'dipoli';  
cfg.tissue = {'brain', 'skull', 'scalp'};  
cfg.conductivity = [0.33 0.0125 0.33];  
headmodel = ft_prepare_headmodel(cfg,  
mesh);
```




Source model computation



```
cfg                = [];  
cfg.elec           = elec_m;  
cfg.headmodel     = headmodel;  
cfg.grid.resolution = 8/1000;  
cfg.grid.unit      = 'm';  
cfg.inwardshift   = 2/1000;  
src_v             = ft_prepare_sourcemodel(cfg);
```

Check if coordinate frames and units are identical

```
figure,  
ft_plot_mesh(headmodel.bnd(1), 'facecolor', 'r',...  
    'edgecolor', 'none', 'facealpha', .1)  
ft_plot_mesh(headmodel.bnd(2), 'facecolor', 'g',...  
    'edgecolor', 'none', 'facealpha', .1)  
ft_plot_mesh(headmodel.bnd(3), 'facecolor', 'b',...  
    'edgecolor', 'none', 'facealpha', .1)  
ft_plot_sens(elec_m, 'label', 'off', 'elecsz', .01,...  
    'elecshape', 'circle', 'facecolor', 'red');  
ft_plot_headshape(hsps, 'vertexcolor', 'm',...  
    'vertexsize', 15)  
ft_plot_mesh(src_v.pos(src_v.inside,:),...  
    'facecolor', 'c', 'vertexsize', 20)  
rotate3d, camlight, view(90,0)
```




Computing the forward solutions

```
cfg                = [];  
cfg.elec           = elec_m;  
cfg.headmodel      = headmodel;  
cfg.grid           = src_v;  
cfg.channel        = raw_clean.label;  
cfg.normalize      = 'yes';  
cfg.backproject    = 'yes';  
cfg.senstype       = 'EEG';  
cfg.unit           = 'm';  
leadfield          =  
ft_prepare_leadfield(cfg, raw_clean);
```

`disp(leadfield)`

```
dim: [17 22 13]  
pos: [4862x3 double]  
unit: 'm'  
inside: [4862x1 logical]  
cfg: [1x1 struct]  
leadfield: {1x4862 cell}  
label: {60x1 cell}  
leadfielddimord: '{pos}_chan_ori'
```



Thanks for your attention!

