

EEG/MEG source reconstruction using beamformers

Amit Jaiswal, DSc (Tech.)

Senior Clinical Scientist & Application Lead

MEGIN Oy, Espoo, Finland



July 25, 2025

IIT Mandi, Mandi, HP, India




MEG/EEG signals measurement


MEG and EEG track electric brain activity by measuring the electromagnetic fields generated by neurons

EEG = measuring the potential differences on the scalp

MEG = measuring neuromagnetic fields outside of the head




MEG + EEG




Source-level analysis: should we do with EEG?

- Due to high spatial resolution MEG source imaging is advantageous in accurately localizing underlying sources.
- Source-level analysis is more common in MEG.
- EEG channels are comparable over participants (within a study) and spatially quite blurry, hence often sensor-level analysis is used.
- EEG source imaging is advantageous if accurate conductor model is provided.




Jaiswal A., 2024

A typical workflow for EEG source imaging




Common source analysis methods

- **Minimum norm estimation (MNE)**
 - estimates the source distribution with the smallest overall power that explains the EEG data.
- **LORETA / sLORETA / eLORETA**
 - methods estimate smooth, distributed source activity.
- **Equivalent current dipole (ECD) fitting**
 - models one or a few discrete dipoles with specific location and orientation.
 - is best for focal sources, such as in epilepsy.
- **Beamforming (e.g., LCMV)**
 - is spatial filtering approach that estimate activity from specific brain regions while suppressing others.
 - Good for analyzing oscillatory activity or task-specific sources.



Beamformers


- Beamforming is scanning method of inverse modeling.
- They scan one source at a time while suppressing activity from other sources (either they are inside brain or noise sources).
- Beamforming is a spatial filtering method, that assumes time courses from difference sources are uncorrelated.
- What it means —
 - The correlated source activity are suppressed.
- Beamformers are adaptive and data-driven, as the filter weights are computed using data covariance.



Sekihara et al., 2001; Spencer et al., 1992; Van Veen et al., 1997, Jaiswal et al., 2020

Common beamformer types

- **LCMV (linearly constrained minimum variance)**
 - Time domain beamformer that minimizes signal power from all directions except the target location.
- **DICS (dynamic imaging of coherent sources)**
 - A frequency-domain beamformer that localizes sources of oscillatory activity.
- **MUSIC (multiple signal classification)**
 - A dipole scanning method that matches observed EEG/MEG signals to theoretical patterns using signal and noise subspaces.




Beamformers in action

Given the sensor-level signal $X(t)$, estimate the activity at source S at location r .

Available ingredient:

- Sensor-level data
- Forward model (lead field)
 - Conductor model
 - Source model
 - Sensor location



Beamformer (LCMV) formulation

Let, x be a $M \times 1$ signal vector of MEG (or EEG) data measured with M sensors, and N be the number of grid points in the source model with grid locations r_j , ($j = 1, \dots, N$). Then the source $y(r_j)$ at any location r_j can be estimated as weighted combination of the measurement x as


$$y(r_j) = W^T(r_j)x$$

where the $M \times 3$ matrix $W(r_j)$ is known as a **spatial filter** for a source at location r_j .

→ *Vector-type beamformer*

Here, $W(r_j)$ is defined as $W(r_j) = \left(L^T(r_j)C^{-1}L(r_j) \right)^{-1} L^T(r_j)C^{-1}$

where $L(r_j)$ is local leadfield matrix, and C is the covariance matrix computed from the measured data samples.



Beamformer formulation

(contd...)

- Since 'noise' typically infects the measured signal, the estimated signal variance is generally normalized using projected noise variance C_n . Such normalized estimates are known as the *Neural Activity Index* (NAI; van Veen et al., 1997) and can be written as

$$NAI(r_j) = \text{Trace} \{ [L^T(r_j)C^{-1}L(r_j)]^{-1} \} / \text{Trace} \{ [L^T(r_j)C_n^{-1}L(r_j)]^{-1} \}$$


- Further, orientation optimization is often used to convert the vector type to a *scalar type beamformer*, turning the spatial filter $W(r_j)$ to $M \times 1$ vector, and computed as

$$w(r_j) = \left(l_{\eta_{opt}}^T(r_j)C^{-1}l_{\eta_{opt}}(r_j) \right)^{-1} l_{\eta_{opt}}^T(r_j)C^{-1}$$

Jaiswal et al., 2020


Best applicability of beamformers

- **Event-related source analysis**
 - Useful for studying task-related brain dynamics by comparing conditions (e.g., rest vs. task).
- **Oscillatory source localization**
 - Ideal for detecting narrow-band neural oscillations (e.g., alpha, beta, gamma) in both MEG and EEG.
- **Exploratory whole-brain scans**
 - Can localize sources without needing prior assumptions about number or location of sources.
- **Good spatial resolution (in MEG)**
 - Particularly effective in MEG due to lower sensitivity to skull conductivity variation compared to EEG.





Limitation of beamformers

- **Poor for highly correlated sources**
 - Beamformers struggle when two or more sources are strongly correlated, often suppressing one or mislocalizing both.
- **Sensitive to head modeling errors**
 - Requires accurate lead field and head model; errors can degrade localization accuracy.
- **Limited for deep sources (especially in EEG)**
 - Deep brain sources are harder to detect, particularly in EEG due to lower sensitivity.
- **Assumes stationarity**
 - Standard beamformers assume that source activity is stationary during the analysis window—limiting for rapidly changing dynamics.



Implementation in open-source s/w and comparison




Jaiswal, A., Nenonen, J., Stenroos, M., Gramfort, A., Dalal, S. S., Westner, B. U., ... & Parkkonen, L. (2020). Comparison of beamformer implementations for MEG source localization. *NeuroImage*, 216, 116797.

Frequency-domain beamformer

- Frequency-domain beamformers localize sources of oscillatory brain activity by applying spatial filters to the cross-spectral density (CSD) of EEG/MEG signals.
- They are ideal for studying frequency-specific activity and functional connectivity in the brain (e.g., alpha, beta rhythms).

Learn more on the method and application on EEG data:

- Groß, J. et al., (2001). Dynamic imaging of coherent sources: studying neural interactions in the human brain. *Proceedings of the National Academy of Sciences*, 98(2), 694-699.
- Tiwari, A. et al. (2022). Dorsal-Ventral Visual Pathways and Object Characteristics: Beamformer Source Analysis of EEG. *Computers, Materials & Continua*, 70(2).



Let's do it in FieldTrip!




Compute data and noise covariance

```
%% Compute noise cov.
```

```
cfg = [];  
cfg.toilim = [-.5, .05];  
epochs_pre = ft_redefinetrial(cfg, epochs);  
cfg = [];  
cfg.covariance = 'yes';  
evoked_pre = ft_timelockanalysis(cfg, epochs_pre);
```

```
%% Compute data cov.
```

```
cfg = [];  
cfg.toilim = [.05, .1];  
epochs_post = ft_redefinetrial(cfg, epochs_pre);  
cfg = [];  
cfg.covariance = 'yes';  
evoked_post = ft_timelockanalysis(cfg, epochs);
```




Computing a common spatial filter

```
cfg = [];  
cfg.method = 'lcmv';  
cfg.grid = leadfield;  
cfg.headmodel = headmodel;  
cfg.lcmv.keepfilter = 'yes';  
cfg.lcmv.fixedori = 'yes';  
cfg.lcmv.reducerank = 3;  
cfg.lcmv.normalize = 'yes';  
cfg.senstype = 'EEG';  
cfg.lcmv.lambda = '5%';  
source_avg =  
ft_sourceanalysis(cfg, evoked);
```

```
disp(source_avg)  
  
time: [1×601 double]  
cfg: [1×1 struct]  
dim: [17 22 13]  
inside: [4862×1 logical]  
pos: [4862×3 double]  
unit: 'm'  
method: 'average'  
avg: [1×1 struct]
```


Apply it noise (baseline) data

```
cfg                = [];  
cfg.method         = 'lcmv';  
cfg.senstype       = 'MEG';  
cfg.grid           = leadfield;  
cfg.grid.filter    = source_avg.avg.filter;  
cfg.headmodel      = headmodel;  
source_pre = ft_sourceanalysis(cfg, evoked_pre); % apply spatial filter on  
noise
```




Apply it to data-of-interest (active time window)

```
cfg                = [];  
cfg.method         = 'lcmv';  
cfg.senstype       = 'MEG';  
cfg.grid           = leadfield;  
cfg.grid.filter    = source_avg.avg.filter;  
cfg.headmodel      = headmodel;  
source_post = ft_sourceanalysis(cfg, evoked_post); % apply spatial filter on  
data
```



Computing NAI (Neural activity index)

```
% Calculate Neural Activity Index (NAI)  
spatial_map = source_post;  
spatial_map.avg.pow =  
(source_post.avg.pow-source_pre.avg.pow)./source_pre.avg.pow;
```



Map the source activity on MRI

```
% Interpolate source map with MRI
```

```
cfg = [];
```

```
cfg.downsample = 2;
```

```
cfg.parameter = 'pow';
```

```
mapOnMRI = ft_sourceinterpolate(cfg, spatial_map , mri);
```

```
% Plot source map on MRI
```

```
maxval = max(mapOnMRI.pow, [], 'all');
```

```
cfg = [];
```

```
cfg.method = 'ortho';
```

```
cfg.funparameter = 'pow';
```

```
cfg.maskparameter = cfg.funparameter;
```


```
cfg.funclorlim = [0.0 maxval];
```

```
cfg.opacitylim = [0.0 maxval];
```

```
cfg.opacitymap = 'rampup';
```

```
cfg.atlas = [ft_dir, 'template/atlas/brainnetome/BNA_MPM_thr25_1.25mm.nii'];
```

```
ft_sourceplot(cfg, mapOnMRI)
```



Thanks for your attention!

